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

The new generation of OPOs

## TUNABLE OPO

Panther<sup>®</sup> EX PLUS OPO  
shown pumped by the Powerlite™ 9010

Panther<sup>®</sup> EX PLUS OPO  
output beam characteristics

#### FEATURES & BENEFITS

*Excellent beam quality and low divergence in both axes*

*Signal energies to > 100 mJ per pulse*

*Linewidths down to < 1.5 cm<sup>-1</sup>*

*Tunability from 205 nm to 2550 nm*

*Pump with any Surelite™ or Powerlite™ pulsed Nd: YAG laser*

*Integrated frequency doubler*

*Hermetically sealed and temperature stabilized crystal housings*

**The Panther<sup>®</sup> EX and Panther<sup>®</sup> EX Plus**, yet another generation of industry leading OPO's.

PANTHER <sup>®</sup> EX & PANTHER <sup>®</sup> EX PLUS SPECIFICATIONS										
	355 nm Laser Output (mJ)	Panther Pump Energy (mJ)	Energy at Peak in Signal unseeded (mJ)	linewidth (cm <sup>-1</sup> )	seeded (mJ)	Energy at Peak in Doubled Signal unseeded (mJ)	linewidth (cm <sup>-1</sup> )	seeded (mJ)	linewidth (cm <sup>-1</sup> )	seeded (mJ)
<b>Panther<sup>®</sup> EX (pumped with 200 mJ or less)</b>										
SL I-10	100	100	12	< 6	N/A	N/A	1.5	< 5	N/A	N/A
SL I-20	100	100	12	< 6	N/A	N/A	1.5	< 5	N/A	N/A
SL II-10	160	160	24	< 6	N/A	N/A	2.5	< 5	N/A	N/A
PL 8030	150	130	20	< 5	25	< 2.5	2.5	< 5	3.5	< 4

Panther <sup>®</sup> EX PLUS (pumped with more than 200 mJ)										
	225	225	35	< 6	N/A	N/A	3.5	< 5	N/A	N/A
SL III-10	225	225	35	< 6	N/A	N/A	3.5	< 5	N/A	N/A
PL 8000	310	250	40	< 5	50	< 2.5	5	< 5	7.5	< 4
PL 9010	550	450	70	< 5	85	< 2.5	7	< 5	10	< 4
PL 9020	475	380	60	< 5	75	< 2.5	6	< 5	9	< 4
PL 9030	400	320	40	< 5	48	< 2.5	5	< 5	7.5	< 4

#### NOTES

All signal specifications at 450 nm and all doubled signal specifications at 240 nm, unless otherwise noted.


As a part of our continuous improvement program, all specifications are subject to change without notice.

## PANTHER® EX OPO & PANTHER® EX PLUS OPO SPECIFICATIONS


DESCRIPTION	SL I,II,III	PL 8000	PL 9000
Pulsewidth (nsec)	3-5	3-7	3-7
Pointing Stability ( $\mu$ rad)	$<\pm 150$	$<\pm 100$	$<\pm 100$
Energy stability (%), 99% of shots)	$<\pm 10\%$ @ 10 & 20 Hz	$<\pm 10\%$ @ 10 & 20 Hz	$<\pm 15\%$ @ 30 Hz
Divergence (mrad, FWHM)		< 2 mrad (both axes)	
Beam Diameter (mm, near field)		4-7	
Beam Ellipticity (near field)		< 20 %	
Polarization			
Signal Vertical %		> 99	
Idler Horizontal %		> 99	

**Continuum®**  
An Excel Technology Company

## PANTHER® EX OPO OPTICAL LAYOUT


## PANTHER® EX PLUS PERFORMANCE WITH SEEDED PL 9010 PUMP

## PANTHER® EX OPO PHYSICAL LAYOUT

All dimensions in mm (inches)

# D Datasheet KIM101<sup>1</sup>

*Appendix C Specifications and Associated Parts*

## Appendix C Specifications and Associated Parts

### C.1 Specifications

Parameter	Value
<b>Piezoelectric Output (SMC)</b>	
Voltage Output	85 to 125 V DC/channel
<b>External Input (SMA CHA, CHB)</b>	
Input Type	Single Ended, Analog
Input Voltage	0 to 10 V ±2% 0 V gives nominally 2kHz pulses backwards 10 V gives nominally 2 kHz pulses forwards 5 V gives zero movement
Output Configuration	Individual Channels: Ch 1, Ch2, Ch3, Ch4 Paired Channels: Ch 1 and 2, Ch 3 and 4
Output Pulses	Frequency: 1 Hz to 2 kHz Voltage: Peak 85 to 125 V Adjustable
Operating Display	128 x 128 TFT LCD Color
I/O 1 and I/O 2 Connectors	TTL Input, TTL Output, 5V Level
Joystick	2-axis, Spring-loaded, Center return Velocity Control of Selected Axis Menu Control in Setup Mode
User I/O Connector	15 Pin DIN 5 V @ 150 mA User Supply with 0 V Return Qty 8 Limit Switch Inputs - Multiplexed
USB Connector	The USB 3.0 (left hand side of the connector is compatible with a USB 2.0 Micro B) Communication speed USB 1.1
<b>Input Power Requirements</b>	
Voltage	15 VDC
Current	2 A
<b>General Data</b>	
Housing Dimensions Excluding Baseplate (W x D x H)	121 x 60 x 47 mm (4.8 x 2.4 x 1.8 in.)
Weight:	390 g

1. Extrait du manuel de KIM101 : [https://www.thorlabs.com/newgroupage9.cfm?objectgroup\\_id=9790&pn=KIM101#9791](https://www.thorlabs.com/newgroupage9.cfm?objectgroup_id=9790&pn=KIM101#9791), consulté le 28/08/2018

# E Datasheet PIAK10<sup>1</sup>

---

## Chapter 5 Specification

Parameter	PIA13	PIA25	PIA50	PIAK10 <sup>a</sup>
Travel Range	13 mm (0.51")	25 mm (0.98")	50 mm (1.97")	10 mm (0.39")
Typical Step Size <sup>b</sup>		20 nm		
Maximum Step Size <sup>c</sup>		<30 nm		
Step Size Adjustability <sup>d</sup>		Up to 30%		
Maximum Step Frequency		2000 Hz		
Velocity (Cont. Stepping)		2 mm/min (Typ), <3.5 mm/min (Max)		
Max Active Preload <sup>e</sup>		25 N	30 N	
Typical Angular Resolution <sup>f</sup>		N/A	1" Mounts: 0.5µrad 2" Mounts: 0.3 µrad	
Recommended Max Active Load Capacity <sup>g</sup>		2.5 Kg		
Mounting	9.5 mm (3/8") Barrel and 3/8-40 Thread		9.5 mm (3/8") Barrel and 1/4-80 Thread	
Drive Screw	1/4-80 Thread, Hard PVD Coated			
Actuator Tip	Tungsten Carbide Ball			
Lifetime	>1 Billion Steps			
Connector Type	SMC Female			
Cable Length	1 m (3.3 ft) Cable Wrap Included			
Cable Exit Adjustability	±55 Deg for Left- or Right-hand Exit			
Piezo Capacitance	175 nF			
Max Operating Voltage	130 V			
Operating Temperature	10 to 40 °C			
Dimensions	59.5 x 31.5 x 17.0 mm (2.34" x 1.24" x 0.67")	71.4 x 31.5 x 17.0 mm (2.81" x 1.24" x 0.67")	96.6 x 31.5 x 17.0 mm (3.80" x 1.24" x 0.67")	72.9 x 31.5 x 17.0 mm (2.87" x 1.24" x 0.67")
Weight (inc. Cable)	55 g (1.94 oz)	60 g (2.12 oz)	65 g (2.29 oz)	55 g (1.94 oz)

Notes.

<sup>a</sup> Kinematic mount version for mirror mounts

<sup>b</sup> Can vary over 20% due to component variance, change of direction and application condition

<sup>c</sup> Adjusted by changing the piezo drive voltage - see the KIM101 controller handbook for more details.

<sup>d</sup> Adjusted by changing the piezo drive voltage - see the KIM101 controller handbook for more details.

<sup>e</sup> Axial force applied onto drive screw tip (ball). A minimum of 5N active preload is recommended to enhance stepping behaviour

<sup>f</sup> Typical angular resolution when fitted to 1" or 2" mirror mounts with max active preload of 30N. Not applicable to PIA13, PIA25 and PIA50 actuators.

<sup>g</sup> A higher horizontal load is possible but this may adversely affect the typical step size.

1. Extrait du manuel de PIAK10 : [https://www.thorlabs.com/newgroupage9.cfm?objectgroup\\_ID=9776](https://www.thorlabs.com/newgroupage9.cfm?objectgroup_ID=9776), consulté le 28/08/2018

# F Datasheet TDS3052S<sup>1</sup>

## Datasheet

### Work where you need to

The TDS3000C Series packs the power of a DPO in a compact design that is only 5.9 in. (149 mm) deep, freeing up valuable benchtop space. And when you need to move your oscilloscope to another lab, its portable 7 lb. (3.2 kg) design makes for easy transport.

If your work demands even more mobility, then the optional battery pack will give you up to three hours of operation without line power.

TDS3BATC provides you with up to three hours of portable battery operation.

## Specifications

All specifications are guaranteed unless noted otherwise. All specifications apply to all models unless noted otherwise.

	TDS3012C	TDS3014C	TDS3032C	TDS3034C	TDS3052C	TDS3054C
Input channels	2	4	2	4	2	4
Bandwidth	100 MHz	100 MHz	300 MHz	300 MHz	500 MHz	500 MHz
Rise time (typical)	3.5 ns	3.5 ns	1.2 ns	1.2 ns	0.7 ns	0.7 ns
Sample rate on each channel	1.25 GS/s	1.25 GS/s	2.5 GS/s	2.5 GS/s	5 GS/s	5 GS/s
Record length	10 kpoints					

### Vertical system

Hardware bandwidth limits	TDS3012C	TDS3014C	TDS3032C	TDS3034C	TDS3052C	TDS3054C
	20 MHz	20 MHz	20 MHz, 150 MHz	20 MHz, 150 MHz	20 MHz, 150 MHz	20 MHz, 150 MHz
Input coupling	AC, DC, GND					
Input impedance	1 MΩ in parallel with 13 pF or 50 Ω					
Input sensitivity range						
1 MΩ	1 mV/div to 10 V/div					
50 Ω	1 mV/div to 1 V/div					
Vertical resolution	9 bits					
Maximum input voltage						
1 MΩ	150 V <sub>RMS</sub> with peaks at ≤400 V					
50 Ω	5 V <sub>RMS</sub> with peaks at ≤30 V					

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**Vertical system**

DC gain accuracy	±2%
Position range	±5 div

**Horizontal system**

Seconds/division range	TDS3012C	TDS3014C	TDS3032C	TDS3034C	TDS3052C	TDS3054C
	4 ns to 10 s	4 ns to 10 s	2 ns to 10 s	2 ns to 10 s	1 ns to 10 s	1 ns to 10 s

Time base accuracy	±20 ppm over any 1 ms time interval
--------------------	-------------------------------------

**Trigger system**

Trigger modes	Auto (supports Roll Mode for 40 ms/div and slower), Normal, Single Sequence
B trigger	Trigger after time or events
Trigger after time range	13.2 ns to 50 s
Trigger after events range	1 to 9,999,999 events
Trigger types	
Edge	Conventional level-driven trigger. Positive or negative slope on any channel. Coupling selections: AC, DC, Noise Reject, HF Reject, LF Reject
Video	Trigger on all lines or individual lines, odd/even or all fields on NTSC, PAL, SECAM
Extended video	Trigger on specific lines in broadcast and non-broadcast (custom) standards and on analog HDTV formats (1080i, 1080p, 720p, 480p). Requires TDS3VID application module
Pulse width (or glitch)	Trigger on a pulse width <, >, =, ≠ to a selectable time limit ranging from 39.6 ns to 50 s
Runt	Trigger on a pulse that crosses one threshold but fails to cross a second threshold before crossing the first again
Rise/fall time	Trigger on pulse edge rates that are either faster or slower than a set rate. Edges can be rising, falling, or either
Pattern	Specifies AND, OR, NAND, NOR when true or false for a specific time
State	Any logic state. Triggerable on rising or falling edge of a clock. Logic triggers can be used on combinations of 2 inputs (not 4)
Comm	Provides isolated pulse triggering required to perform DS1/DS3 telecommunications mask testing per ANSI T1.102 standard. Requires TDS3TMT application module
Alternate	Sequentially uses each active channel as a trigger source

**Acquisition system**

DPO	Captures and displays complex waveforms, random events and subtle patterns in actual signal behavior. DPOs provide 3 dimensions of signal information in real time: Amplitude, time, and the distribution of amplitude over time
Sample	Sample data only
Average	Waveform averaged, selectable from 2 to 512
Envelope	Min-max values acquired over one or more acquisitions
Peak detect	High-frequency and random glitch capture. Captures glitches as narrow as 1 ns (typical) using acquisition hardware at all time base settings

## Datasheet

### Acquisition system

WaveAlert®	Monitors the incoming signals on all channels and alerts the user to any waveform that deviates from the normal waveform being acquired
Single sequence	Use the Single Sequence button to capture a single triggered acquisition sequence at a time

### Waveform measurements

Cursors	Amplitude, time
Automatic measurements	Display any four measurements from any combination of waveforms. Or display all measurements with measurement snapshot feature. Measurements include Period, Frequency, +Width, -Width, Rise time, Fall time, +Duty cycle, -Duty cycle, +Overshoot, High, Low, Max, Min, Peak-to-peak, Amplitude, Mean, Cycle mean, RMS, Cycle RMS, Burst width, Delay, Phase, Area <sup>1</sup> , Cycle Area <sup>1</sup>
Measurement statistics	Mean, Min, Max, Standard deviation. Requires TDS3AAM application module
Thresholds	User-definable thresholds for automatic measurements; settable in percent or voltage
Gating	Isolate a specific occurrence within an acquisition to take measurements, using either the screen or cursors

### Waveform math

Arithmetic	Add, subtract, multiply, and divide waveforms
FFT	Spectral magnitude. Set FFT vertical scale to Linear RMS or dBV RMS, and FFT window to Rectangular, Hamming, Hanning, or Blackman-Harris
Advanced math	Integrate, differentiate, define extensive algebraic expressions including analog waveforms, math functions, scalars, up to two user-adjustable variables and results of parametric measurements. For example: (Intg (Ch1-Mean(Ch1)) × 1.414 × VAR1) <sup>2</sup>

### Waveform processing

Autoset	Single-button, automatic setup of all channels for vertical, horizontal and trigger systems, with undo autoset
Deskew	Channel-to-channel deskew ±10 ns may be manually entered for better timing measurements and more accurate math waveforms

### Display system

Display type	6.5 in. (165.1 mm) liquid-crystal TFT color display
Display resolution	640 horizontal × 480 vertical pixels (VGA)
Interpolation	Sin(x)/x
Waveform styles	Dots, vectors, variable persistence, infinite persistence
Graticules	Full, grid, crosshair, and frame. NTSC, PAL, SECAM, and vectorscope (100% and 75% color bars) with optional TDS3VID application module
Format	YT, XY, and Gated XYZ (XY with Z-axis blanking available on 4-channel instruments only)

<sup>1</sup> Requires TDS3AAM application module.

<sup>2</sup> Requires TDS3AAM application module.

**Input/output ports**

<b>Ethernet port</b>	RJ-45 connector, supports 10Base-T LAN
<b>USB port</b>	Front-panel USB 2.0 host port. Supports USB flash drive
<b>GPIB port</b>	Full talk/listen modes, setting and measurements. (Optional with TDS3GV Communications Module)
<b>RS-232C port</b>	DB-9 male connector, full talk/listen modes; control of all modes, settings and measurements Baud rates up to 38,400 (Optional with TDS3GV Communications Module)
<b>VGA video</b>	DB-15 female connector, monitor output for direct display on large VGA-equipped monitors (Optional with TDS3GV communications module)
<b>External trigger input</b>	BNC connector, input impedance >1 MΩ in parallel with 17 pF; max input voltage is 150 V <sub>RMS</sub>

**Power source**

<b>AC line power</b>	
<b>Source voltage</b>	100 V <sub>RMS</sub> to 240 V <sub>RMS</sub> ±10%
<b>Source frequency</b>	45 Hz to 440 Hz from 100 V to 120 V 45 Hz to 66 Hz from 120 V to 240 V
<b>Power consumption</b>	75 W maximum
<b>Battery power</b>	Requires TDS3BATC, rechargeable lithium ion battery pack
<b>Operating time, typical</b>	3 hours

**Physical characteristics**

Dimensions		mm	inches
	Width	375.0	14.8
	Height	176.0	6.9
	Depth	149.0	5.9
Weight		kg	lb.
	Instrument only	3.2	7.0
	with accessories	4.5	9.8
Package dimensions		mm	inches
	Width	502.0	19.8
	Height	375.0	14.8
	Depth	369.0	14.5
Rackmount (RM3000)		mm	inches
	Width	484.0	19.0
	Height	178.0	7.0
	Depth	152.0	6.0

# G Datasheet DET10A <sup>1</sup>

Si Biased Detector

Chapter 7: Specifications

## Chapter 7 Specifications

All measurements are performed at 25 °C ambient temperature unless stated otherwise.

Electrical Specifications		
<b>Detector</b>	-	Silicon
<b>Active Area</b>	-	Ø1.0 mm (0.8 mm <sup>2</sup> )
<b>Wavelength Range<sup>a</sup></b>	$\lambda$	200 to 1100 nm
<b>Peak Wavelength</b>	$\lambda_p$	730 nm (Typ.)
<b>Peak Response<sup>b</sup></b>	$R(\lambda_p)$	0.44 A/W (Typ.)
<b>Shunt Resistance</b>	$R_{sh}$	1 GΩ (Typ.)
<b>Diode Capacitance</b>	$C_J$	6 pF (Typ.)
<b>Rise Time<sup>b,c,d</sup></b>	$t_r$	1 ns (Typ.)
<b>NEP (830 nm)</b>	-	$5.0 \times 10^{-14} \text{ W}/\sqrt{\text{Hz}}$ (Typ.)
<b>Photodiode Bias Voltage</b>	$V_R$	10 V
<b>Dark Current<sup>e</sup></b>	$I_D$	0.3 nA (Typ.) 2.5 nA (Max)
<b>Output Voltage</b>	$V_{OUT}$	0 to 10 V
General		
<b>On/Off Switch</b>	Slide	
<b>Battery Check Switch</b>	Momentary Pushbutton	
<b>Output</b>	BNC (DC Coupled)	
<b>Package Size</b>	2.8" x 1.9" x 0.83" (70 mm x 48 mm x 21 mm)	
<b>PD Surface Depth</b>	0.13" (3.4 mm)	
<b>Weight</b>	0.2 kg	
<b>Accessories</b>	SM1T1 Coupler SM1RR Retainer Ring	
<b>Storage Temp</b>	-20 to 70 °C	
<b>Operating Temp</b>	10 to 50 °C	
<b>Battery</b>	A23, 12 V <sub>DC</sub> , 40 mAh	
<b>Low Battery Voltage<sup>f</sup></b>	(See Battery Check)	
$V_{OUT}$ (Hi-Z)	~9 V	
$V_{OUT}$ (50Ω)	~400 mV	

<sup>a</sup> When long-term UV light is applied, the product specifications may degrade. For example, the product's UV response may decrease and the dark current may increase. The degree to which the specifications may degrade is based upon factors such as the irradiation level, intensity, and usage time.

<sup>b</sup> Measured with specified bias voltage of 10.0 V.

<sup>c</sup> Low battery voltage will result in slower rise times and decreased bandwidth.

<sup>d</sup> For a 50 Ω Load

<sup>e</sup> Measured with a 1 MΩ Load

<sup>f</sup> Assumes the battery voltage drops below 9.6 V. The reverse protection diode generates a 0.6 V drop.

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<sup>1</sup> For higher loads, the flip time must be increased. Settings will differ for each application, but as a rough guide, for loads >50 g the flip time should be > 1 second.

1. Extrait du manuel de MFF101R : [https://www.thorlabs.com/newgroupage9.cfm?objectgroup\\_id=3962&pn=MFF101/M#3962](https://www.thorlabs.com/newgroupage9.cfm?objectgroup_id=3962&pn=MFF101/M#3962), consulté le 28/08/2018

# I Datasheet PM1918R<sup>1</sup>

---

22

General Information

and manufacturing environments, the 1918-R Series optical meters enable users to benefit from high correlation between multiple locations at a price-to-performance ratio second to none. Among all the other practical tools provided, remote controlling with a computer and synchronization to other instruments are simplified with the inclusion of LabVIEW drivers, tools to develop in the .NET environment and a high-speed software utility that fully utilized the optical meter's ability to sample at 250 kHz and transfer data via a Full-Speed USB interface (11 Mbps).

## 2.2 Calibration

---

Calibration of the optical meter is done at the factory by defining a slope and offset, for all ranges as determined for each detector type.

Newport recommends annual factory re-calibration to ensure the continued accuracy of optical meter measurements.

Please refer to the "Maintenance and Troubleshooting" section for contact information for re-calibration of your optical meter.

## 2.3 Specifications

---

### 2.3.1 Physical Specifications:

Dimensions:	5.4 x 7.6 x 2.4 inch (137 x 193 x 61 mm)
Weight:	max. 2.3 lb
Enclosure:	Plastic molded case with a kick stand.
Connectors:	Optical Detector 15-Pin D-Sub, Analog Output 3.5 mm Jack, USB-A Host, mini-B USB Device, DC Input 2.1 mm Power Jack.
Power:	90-264 VAC, 50/60 Hz, 33 Watts with 1918-EXT-PS External Power Supply
Display:	Graphical LCD ¼ VGA, 4 inch diagonal
Display Update Rate:	up to 20 Hz for selected graphical modes
Operating Environment <sup>1</sup> :	5°C to 40°C; < 70% RH non-condensing
Storage Environment:	-20°C to 60°C; < 90% RH non-condensing
Altitude	< 3000m
Installation Category	II
Pollution Degree	2
Use Location	Indoor use only.

---

<sup>1</sup> Charging battery up to 35°C

### 2.3.2 Electrical Specifications

#### DC Current Measurement (Low-Power, Semiconductor Photodiode)

Signal Range <sup>2,3</sup>	0	1	2	3	4
Full-Scale Current	2.50 $\mu$ A	25.0 $\mu$ A	250 $\mu$ A	2.50 mA	25 mA
Resolution (Filtered) <sup>1</sup>	10 pA	100 pA	1.00 nA	10.0 nA	100 nA
Bandwidth <sup>6</sup> (Unfiltered)	480 Hz	4 kHz	43 kHz	170 kHz	500 kHz
Accuracy (Filtered) <sup>1</sup>	0.2%	0.2%	0.2%	0.2%	0.2%
Accuracy (Maximum Measurement Rate=4 kHz)	1.0%	1.0%	1.0%	1.0%	1.0%
Signal to Noise Ratio <sup>1</sup>	100 dB	100 dB	100 dB	100 dB	100 dB

#### Peak-Peak Current Measurement (Low-Power, Semiconductor Photodiode)

Signal Range <sup>2,3</sup>	0	1	2	3	4
Full-Scale Current	2.50 $\mu$ A	25.0 $\mu$ A	250 $\mu$ A	2.50 mA	25 mA
Resolution	76.3 pA	763 pA	7.63 nA	76.3 nA	763 nA
Accuracy	1.0%	1.0%	1.0%	1.0%	1.0%
Bandwidth <sup>6</sup>	480 Hz	4 kHz	43 kHz	170 kHz	500 kHz
Maximum Pulse Repetition Rate <sup>4</sup>	480 Hz	4 kHz	4 kHz	4 kHz	4 kHz

#### D.C. Voltage Measurement (Thermopile)

Signal Range <sup>2,3</sup>	0	1	2	3	4	5
Full-Scale Voltage	2.5mV	25.0 mV	250 mV	2.50 V	25.0 V	130 V
Resolution	76.3 nV	763 nV	7.63 $\mu$ V	76.3 $\mu$ V	763 $\mu$ V	3.96 mV
Accuracy (Filtered) <sup>1</sup>	0.2%	0.2%	0.2%	0.2%	0.2%	0.2%
Accuracy (Maximum Measurement Rate=2 kHz)	1.0%	1.0%	1.0%	1.0%	1.0%	1.0%
Bandwidth <sup>6</sup>	3 kHz	3 kHz	3 kHz	3 kHz	3 kHz	3 kHz
Signal to Noise Ratio <sup>1</sup>	100dB	100dB	100dB	100dB	100dB	100dB

#### Pulse Voltage Measurement (Energy, Pyroelectric)

Signal Range <sup>2,3</sup>	0	1	2	3	4	5
Full-Scale Current	2.5mV	25.0 mV	250 mV	2.50 V	25.0 V	130 V
Resolution	76.3 nV	763 nV	7.63 $\mu$ V	76.3 $\mu$ V	763 $\mu$ V	3.96 mV
Accuracy	1.0%	1.0%	1.0%	1.0%	1.0%	1.0%
Bandwidth <sup>6</sup>	600 kHz	1000 kHz	1000 kHz	1000 kHz	400 kHz	1000 kHz
Maximum Pulse Repetition Rate <sup>5</sup>	2 kHz	2 kHz	2 kHz	2 kHz	2 kHz	2 kHz

**RMS Measurement**

Maximum Pulse Repetition Rate	2 kHz	4 kHz
Accuracy	1.0%	2.0 %

**Analog Output**

Output Range	0	1	2
Full Scale Voltage (Load > 100 kΩ)	1 V	2 V	5 V
Accuracy	1.0 %	1.0 %	1.0 %
Linearity	0.3%	0.3%	0.3%

**Trigger Level**

Programmable Level	0...100 % Full Scale
Resolution	0.39 % Full Scale

<sup>1</sup> With 0.5 HZ filter on.<sup>2</sup> Listed signal ranges specify meter capability. Available signal ranges are detector dependent.<sup>3</sup> Maximum measurable signal is detector dependent.<sup>4</sup> While the maximum repetition range may equal the bandwidth, it really depends on the signal duty-cycle or the signal shape.<sup>5</sup> The Maximum Repetition Rate refers to the meter pulse-by-pulse measuring capability. Due to its high bandwidth, the 1918-R can take in signals with higher repetitive rates and outputs them undistorted at the Analog Output.<sup>6</sup> The instrument bandwidth is determined by the detector used. Please refer to Newport Corporation's complete offering on detector type. The specified bandwidth is measured from the instrument input (detector) to the Analog Output BNC, using a generator with 50 Ω output.

---

**2.4 Unpacking and Handling**

It is recommended that the 1918-R Optical Meter be unpacked in a lab environment or work site. Unpack the system carefully; small parts and cables are included with the instrument. Inspect the box carefully for loose parts before disposing of the packaging. You are urged to save the packaging material in case you need to ship your equipment in the future.

---

**2.5 Inspection for Damage**

The 1918-R Optical Meter is carefully packaged at the factory to minimize the possibility of damage during shipping. Inspect the box for external signs of damage or mishandling. Inspect the contents for damage. If there is visible damage to the instrument or accessories upon receipt, inform the shipping company and Newport Corporation immediately.

# J Datasheet KDC101<sup>1</sup>

---

## Appendix C Specifications and Associated Products

### C.1 Specifications

Parameter	Value
<b>Motor Output</b>	
Motor Drive Voltage	±12 to ±15V (Depending on Supply)
Motor Drive Current	150mA (Cont) >250 mA (peak)
Motor Drive Type	8-bit Sign/Magnitude PWM
Control Algorithm	Digital PID Filter (16bit)
Position Feedback: Quadrature Encoder (QEP) Input	5V Single Ended
Encoder Feedback Bandwidth	750 kHz
Position Counter	32-bit
Operating Modes	Position, Velocity
Velocity Profile	Trapezoidal
<b>Motor Drive Connector (15 Way D-Type)</b>	
Motor Drive Outputs	+ve & -ve
Quadrature Encoder (QEP) Input	Single Ended
Limit Switch Inputs	Forward, Reverse (+ Common Return)
Encoder Supply	5V
<b>Front Panel Controls</b>	
Sprung Potentiometer Wheel	Bidirectional Velocity Control, Forward/Reverse Jogging or Position Presets
<b>Input Power Requirements</b>	
Voltage	15V Regulated DC
Current	500mA (peak)
<b>General Data</b>	
Housing Dimensions (W x D x H) (excluding buttons and baseplate)	60 x 60 x 47mm (2.36" x 2.36" x 1.85")
Weight	170g (6.0 oz)

---

1. Extrait du manuel de KDC101 : [https://www.thorlabs.com/newgroupage9.cfm?objectgroup\\_id=2419&pn=KDC101#5077](https://www.thorlabs.com/newgroupage9.cfm?objectgroup_id=2419&pn=KDC101#5077), consulté le 28/08/2018

*Appendix C*

---

Recommended Motor Requirements	
Peak Power	2.5W
Rated Current	10mA to 200mA (Nominal)
Motor Type	Brushed DC
Coil Resistance	5 to 50Ω
Coil Inductance	250 to 1500mH
Position Control	Closed loop Encoder

## C.2 Associated Products

Product Name	Part Number
6mm DC Servo Motor Actuator, 1/4"-80	Z806
6mm DC Servo Motor Actuator, 1/4"-80, Vacuum Rated	Z806V
12mm DC Servo Motor Actuator, 1/4"-80,	Z812
12mm DC Servo Motor Actuator, 3/8" Barrel Attachment	Z812B
12mm DC Servo Motor Actuator, 1/4"-80, Vacuum Rated	Z812V
12mm DC Servo Motor Actuator, Vacuum Rated, 3/8" Barrel Attachment	Z812BV
25mm DC Servo Motor Actuator, 3/8" Barrel Attachment,	Z825B
25mm DC Servo Motor Actuator, Vacuum Rated, 3/8" Barrel Attachment	Z825BV
3-Channel K-Cube Controller USB Hub	KCH301
6-Channel K-Cube Controller USB Hub	KCH601
Single Way Power Supply	KPS101
8-way Power Supply	TPS008

# K Datasheet PRMTZ8M<sup>1</sup>

---

## Chapter 5 Specifications

Parameter	Value
<b>Stage Specifications</b>	
Range	360° Continuous
Bidirectional Repeatability	±0.1°
Backlash	±0.3°
Max Rotation Velocity	25 deg/sec
Calculated Resolution	2 arcsec (0.0005°)
On-Axis Load Capacity, Horizontal Mounting <sup>a</sup>	15 N (1.5 kg/3.4 lbs)
Min Repeatable Incremental Motion	0.04°
Absolute On-Axis Accuracy	0.1%
Home Location Accuracy	±0.2°
Wobble	200 µrad
<b>Motor Specifications</b>	
Motor Type	DC Servo
Nominal voltage	6 V
No load speed	6560 rpm
No load current	7.43 mA
Nominal speed	1050 rpm
Nominal torque (max. continuous torque)	1.26 mNm
Nominal current (max. continuous current)	0.156 A
Stall torque	1.54 mNm
Starting current	0.184 A
Max. efficiency	65%
Terminal resistance	32.7 Ω
Terminal inductance	0.607 mH
Torque constant	8.38 mNm/A
Speed constant	1140 rpm/V
Speed / torque gradient	4450 rpm/mNm
Mechanical time constant	13.9 ms
Rotor inertia	0.298gcm <sup>2</sup>

<sup>a</sup> The stage is not intended for vertical mounting, e.g. on a post.

1. Extrait du manuel de PRMTZ8M : [https://www.thorlabs.com/newgroupage9.cfm?objectgroup\\_id=4134&pn=PRMTZ8/M](https://www.thorlabs.com/newgroupage9.cfm?objectgroup_id=4134&pn=PRMTZ8/M), consulté le 28/08/2018

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